ProMotion®

INNOVATIVE
THERAPY ON
THE MOVE

ProMetion

Official Lead Sponsor

FOOTBALL MEDICINE & PERFORMANCE ASSOCIATION

WWW.PROMOTION.FITNESS



CONTENTS



ABOUT US



PROMOTION EV1



FOR ATHLETES



FOR PHYSIOTHERAPISTS

THE SCIENCE



ELITE PACKAGE



ProMOTION EV1 is cutting edge technology created by Swellaway Limited to support athletes to get back to fitness, faster. The Swellaway team has worked closely with elite sports clubs, physiotherapists and University academics to develop this innovative new technology.

Each stage of product development has been backed by robust scientific evidence thanks to a unique Knowledge Transfer Partnership (KTP) with University of Central Lancashire and Manchester Metropolitan University.







ProMoTION

INNOVATIVE THERAPY ON THE MOVE

Through these partnerships, we've worked hard to truly understand the therapy challenges that athletes and physiotherapists face.

ProMOTION EV1 helps to solve these challenges:



REDUCE RECOVERY TIME AND GET BACK TO FITNESS FASTER



CORDLESS PORTABLE DEVICE FOR USE ON THE MOVE



TOOLS TO MONITOR USAGE OUTSIDE THE TREATMENT ROOM



PRECISE TEMPERATURE AND COMPRESSION CONTROL



TARGETED THERAPY TO MINIMISE REDUCTION IN MUSCLE STRENGTH

ProMoTION[®]

PROMOTION EVII

O PORTABLE

Portable device to allow for therapy on the move

Every aspect of ProMOTION EV1 is designed to get athletes back to fitness faster.



Pre-set programs for players to use at home



History logs to track therapy usage



To minimise reduction in muscle strength



Temperature range from 6-40°c, with rapid temperature changes

OMPRESSION RANGE

Compression range from 20mmHg-75mmHg

PRECISE CONTROL

Precise maintenance of temperature and compression

MANAGE THERAPY

Manage therapy by creating your own pre-set programs





FOR PHYSIOTHERAPISTS

COMPLEMENT YOUR EXISTING TOOLS

ProMOTION EV1 bolsters the range of therapy options you can offer to support athletes back to fitness.

This unique wellness device provides hot, cold or contrast therapy, between the ranges of 6°C – 40°C and the option of compression from 20mmHg to 75mmHg. It's fully portable with no need for ice, water or tethering.

"ProMOTION EV1 is an exciting product for therapists. To have the ability to treat patients passively with a hand-held unit, whilst also having the option to allow them to apply the same unit without the need for a therapist is unique to the market."

CARL NELSON

PREMIER LEAGUE PHYSIOTHERAPIST

"I WOULD DEFINITELY HAVE A PROMOTION EV1 IN MY TOOLKIT."

NICK MURPHY

WARRINGTON WOLVES

"ProMOTION EV1 allows us as a medical team to focus on particular areas of the body to provide controlled temperature and compression. ProMOTION EV1 is a great device to help our players be refreshed and fully recovered in the shortest possible time."

JACK TARGETT

BRISTOL BEARS RUGBY



PERSONALISED PROGRAMS

Create personalised programs with specific temperature levels, compression levels and duration.

MANAGE THERAPY

Manage therapy by creating your own pre-set programs.

HISTORY LOG

Create a personal or organisational profile to log your injury or recovery information.

SHARE HISTORY

Share your history with others.







INNOVATION BACKED BY SCIENCE

ProMOTION EV1 is the first product based on new technology and science developed alongside academics from the University of Central Lancashire and Manchester Metropolitan University.

Through our Knowledge Transfer Partnership, the device has been through years of trials and rigorous testing with professional teams and leading academics to optimise safety and performance.

WHY TARGETED CRYOTHERAPY?

The guidelines for injury management have introduced 'optimal loading' as a solution to encourage early activity to optimise recovery.

Evidence supports the use of cryotherapy as short-term pain relief for acute soft-tissue injury management. However, short-term adverse effects to dynamic stability and muscle strength (up to a 16% reduction) have also been reported following global ice interventions. Therefore, authors have warned clinicians to consider the potential increased risk of injury following a global cryotherapy treatment, before returning to weight bearing activity.

Emerging research has now shown that targeted cryotherapy, such as that offered by *ProMOTION* EV1, has the potential to provide necessary relief and minimise adverse effects. Athletes can then return to weight-bearing activity earlier.



THE FLORE PACKAGE



WHAT'S INCLUDED:



EVI THERAPY UNIT



COMPRESSION PUMP



BATTERIES



COMPRESSION WRAPS



UNIVERSAL STRAP



PROMOTION APP



BATTERY CHARGER



CARRY CASE

FOR MORE INFORMATION

GET IN TOUCH

Jonathan Smith

Managing Director 07990 586 290 jsmith@swellaway.com

Rebecca Woodley-Jones

Commercial Director 07715 036 546 rebecca@swellaway.com

Olivia Greenhalgh

Sports Rehabilitation Technologist 01925 751 027 olivia@swellaway.com

General Enquiries

For all other enquiries 01925 751 027 customerservices@swellaway.com



